

Technical Specifications (Technical specifications subject to change without notice!)

Specifications	Model FIXI	Model FOX
Weight / product (kg)	6.2	6.3
Max. body weight (kg)	100	100
Length approx.	680	680
Width approx.	580	580
Seat height (mm)	470	465
Seat width (mm)	380	385
handle height (mm)	610-700 / 5 positions	700-780 / 5 positions
Turning radius approx. (mm)	900	900
wheels	4x Plain bearing, D=200	4x Plain bearing, D=200
Type of brake	2x permanent brakes	2x permanent brakes
Angle between longitudinal centre line of handgrip and movement direction: 0 degrees (parallel)		

* Weight only accessories Production-dependent tolerances of tubes +/- 15 mm

Spare Parts and Accessories

Art. No.	name	Description
288.00.45	Front wheel, black	plastic
288.00.46	Rear wheel, black	plastic
280.00.60	Brake set 1.265x1.000 Niro incl. nipple	VA steel
280.10.48	Black cushion (accessories: weight 70 g)	PU
286.65.25	Seat, black (accessories: weight 800 g)	wood

Advisory!

Always keep the User Instructions in case of future use.

Always bring a copy of the User Instructions with you whenever you leave the house so that in the case of an emergency another person will be able to operate the Rollator.

The identical serial number in the User Instructions is also printed on the nameplate of the product. For complaints or questions, always provide the model, serial no. and the User Instructions no.

The user instructions and the maintenance manual are downloadable on our homepage.



User Instructions

BA No.: 02.05.012/2016-12



REBOTEC

REBOTEC® Rehabilitationsmittel GmbH

PO Box Q541
Queen Victoria Building
NSW 1230

rebotec.com.au
info@rebotec.com.au



Rollator For Kids

FIXI
186.00.10
Serial #

FOX
186.00.20
Serial #



With seat and cushion (accessories)

Intended Use / User Advisory

The Rollators are designed for persons with limited walking ability and serve to increase their ability to move. The Rollators may be used only on solid, even surfaces.

The aluminium tubular design provides security against tipping over. Its lightweight not only makes the Rollator easier to drive, the bow handle also allows it to be pulled up towards the upper body, thus allowing the patient to easily overcome street kerbs. In addition the backrest may be outfitted with a back padding (accessory), which makes sitting comfortable. For transport, the Rollator can be folded up. The bow handle is also height adjustable to suit the height of the patient.

Have your specialised dealer or care provider advice you on using the Rollator and its adjustment for proper function!

Instructions for Assembly and First-time Use

1. Unfolding the Rollator

Place the Rollator on its front wheels and push the frame with the rear wheels apart (A) **Fig. 1**. The Rollator clicks by itself into its final position (B) **Fig. 2**.

2. Adjusting the Handles-Fig. 3

Next remove both star-grip screws (C). Then choose the appropriate height (D) and adjust the grip handles to the next screw hole. Reinstall the star-grip screws (C) and tighten them firmly.

3. Using the Brakes –Fig. 4

By pulling up the left and right brake levers (G), the back wheels are locked respectively. The permanent brake is set by pushing the brake levers (H) downwards. By pulling up the brake levers (G) the permanent brake is disengaged.

4. Folding the Rollator together

Then with one hand pull the **red cord** (J) **Fig. 2** upwards until the click-lock disengages. Then with the other hand lift the Rollator at the tubular bow (F).

5. Using the Seat (accessory)-Fig. 2

The seat (M) is screwed firmly to the frame (F).



- Before using the seat, activate the permanent brake (H), see Step 4.
- When the sun shines the seat may get hot! Therefore, check the temperature of the seat with your hand. If necessary, allow the seat to cool down, or place a cushion on the seat.

6. Adjusting the Brakes-Fig. 5

Loosen the counter nuts (L) and then adjust the brake with the adjustment screw (K). Be sure that when the brake is released, the wheels are allowed to turn freely, and when the brake levers are pulled the wheels are firmly arrested. Then retighten the counter nut (L).



Safety Instructions and Inspection

- After the first time set-up and before each use, always conduct a visual inspection of the Rollator for deficiencies. Check the screw connections to be sure they are tight. Tighten if necessary using the appropriate tools. If you have any doubt, then do not use the Rollator. Contact a specialised dealer.
 - Before each use, check over the folding mechanism to be sure it has clicked into its final position.
 - Check the function of the brakes and the permanent brake function as well.
 - Use the grasp handle only for raising the Rollator and as a backrest.
 - **DANGER OF TIPPING!**
- For moving the Rollator only use the rubberized handles.

Care and Maintenance

The surface of the aluminium tubing has been provided with a paint coating. The Rollator can be cleaned with a commercially available cleaning or disinfecting agent – always follow instructions for use. All moving parts are permanently lubricated. **Inadequate or negligent care and maintenance leads to limitations or annulment of product liability by the Manufacturer.**

Allow your Rollator to be checked over by specialised dealer once a year!

Reuse

The product is intended for reuse. The reuse (life) depends on the product state. Before the re-use of the device maintenance must be carried out and the appropriate instructions to security, care and maintenance must be observed. This review shall exclude a possible improper use caused by obvious or invisible damage. Please refer to the Medical Device Directive (MDD) and the regulation for operator (MPBetreibV).

Life

The expected life of the products listed in this manual when used in accordance with safety instructions, maintenance intervals and care instructions is 3 years. Beyond this period, the product can be set used only in a safe condition. Maintenance should be done at least every six months. Specific requirements of the product, due to the use of the patient, causing shorter maintenance intervals.

Warranty

The products listed in this manual are provided with a warranty of 24 months. Wear parts and damages of the products due to improper use are excluded.

Figures:

