

**Important Note:** These instructions are general information only and do not take into account your specific circumstance. These instructions are to be used in conjunction with the information provided by your health professional. Do not exceed the Maximum User Weight specified by the manufacturer. If you are unsure on how to use or fit a piece of equipment, please contact your Health Professional.

## Safety Information and Precautions

- Dual motor lift recline chairs are designed to provide comfort and ease of movement for users, offering independent control of the backrest and footrest for optimal positioning. It is important to understand the correct use and care of the chair to ensure safety and maintain its condition.



### Using a Lift Recline Chair

- Before use, ensure the chair properly positioned, with sufficient room to operate without coming into contact with walls or other objects.
- Ensure it is plugged in and all controls are functioning properly.
- Use the remote control to adjust the backrest and footrest independently. Follow the user manual to familiarise yourself with the correct operation of the chair.
- Ensure that the floor area around the chair is clear of obstacles to avoid tripping hazards during transfers or when adjusting the chair.
- Keep the chair on a level surface to ensure stability when raising or reclining.
- If the chair has wheels, ensure brakes are applied.

### Precautions and Safety

- Do not exceed the Safe Working Load as specified by the manufacturer.
- Always ensure the user is fully seated and stable before operating the recline or lift functions. Be cautious when using the lift function to avoid tipping or sudden movements that may cause injury.
- Do not allow children to operate the chair's controls or play around the moving parts.
- Check for pets under or around the chair before lowering or raising the chair to prevent injury or damage.
- Be mindful of pets scratching or biting the leather or fabric upholstery, which can cause permanent damage to the chair's surface.
- Ensure hands, arms, and legs are kept clear of moving parts when adjusting the chair to avoid injury.

- Regularly inspect the power cord and electrical components for any damage. If damage is found, discontinue use and contact us immediately.
- Do not force the chair to move beyond its designed range of motion. Always allow the motors to complete their actions before adjusting further.

### Maintenance and Cleaning

- Clean the leather or fabric regularly according to the manufacturer's instructions to maintain its appearance and prevent wear.
- If the chair is upholstered in leather, use a leather conditioner to prevent cracking or drying out over time.
- Protect the chair from sharp objects such as pet claws, keys, or belts, which can cause tears or scratches in the upholstery.
- Keep liquids and food away from the chair to avoid spills that may damage the internal mechanisms or stain the upholstery.
- Regularly check the chair for loose bolts, wear on electrical components, or other signs of damage. If any issues are found, do not use the chair until it has been repaired.

**Disclaimer:** This product safety information sheet is provided for general guidance only and does not replace the manufacturer's user manual or specific safety protocols recommended by the manufacturer. Always adhere to local safety regulations and the manufacturer's recommendations for maintenance, inspection, and operation. Failure to follow proper procedures and obtain appropriate training may result in injury or damage to the equipment.