

# DELUXE

## - Aluminium Wheelchair


- Luxury pack - semi-adaptable to suit the patient.
- Suitable for hospitals, nursing homes and other institutional use
- Equipped with an array of important features.
- Heavy duty construction.
- Tyres crafted from durable, soft polyurethane.
- No tyre inflation required.
- Designed to fit many optional accessories
- Suitable for domestic use.



Available with either a seat width 46cm or 51cm



Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166  
 22 / 170 North Road, Woodridge, QLD 4119 | **Ph:1300 460 070**  
 wecare@mobilityrentals.com.au | www.mobilityrentals.com.au



### seat width 51cm

overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	backrest height	footrest length	folded width	rear wheel	front castor	unit weight	weight capacity
min90cm max98cm	107cm	68cm	44cm	51cm	51cm	24cm	43cm	min30cm max40cm	30cm	propel 60cm transit 30cm	18cm	16.9kg	136kg

### seat width 46cm

overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	backrest height	footrest length	folded width	rear wheel	front castor	unit weight	weight capacity
min90cm max98cm	107cm	63cm	44cm	46cm	51cm	24cm	43cm	min30cm max40cm	30cm	propel 60cm transit 30cm	18cm	16.9kg	136kg