

# PREMIUM PRESCRIPTION WHEELCHAIR - Aluminium



- Premium pack - Fully adaptable prescription wheelchair.
- Designed to be multi-adjustable.
- Can be tailored for the individual and with a host of accessory choices.
- Equipped with an array of important features as standard.
- Heavy duty construction.
- Available with a choice of rims and tyres.

Available seat widths  
**40cm or 46cm or 51cm**

- (M.) Handle bars**
  - Height adjustable handles
- (E.) Backrest**
  - tension straps
  - Height adjustable
  - Thick padded cover
  - Spill proof cloth
  - hypoallergenic
- (F.) Backrest**
  - angle adjustable
- (G.) Armrest**
  - flip back motion
  - removable
  - height adjustable
  - slide mechanism armpad (forward/backward)
  - padded
- (B.) Seat**
  - Depth adjustable
- Seat**
  - tension adjustable straps
  - standard cushion included
  - Spill proof cloth (black)
  - hypoallergenic
- (H.) Legrest**
  - swings inward and outward
  - height adjustable
- (I.) Footplate**
  - height adjustable
  - angle adjustable
  - fore and aft position adjustable
  - heel straps
- (L.) Castors**
  - standard castors or
  - suspension castors as optional
- (D.) Castors**
  - adjustable seat height
- (K.) Tyres**
  - pneumatic tyres as standard or
  - "soft ride" polyurethane as optional
- Spoked wheels**
- Quick release wheels**
- Park brake**

Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166  
 P O Box 249 Underwood QLD 4119 | **Ph: 1300 460 070**  
 wecare@mobilityrentals.com.au | www.mobilityrentals.com.au



**1 YEAR WARRANTY**  
\* conditions apply

- (D.) Axle plate** to adjust seat height
- (J.) Axle plate** to adjust wheel base
- Anti-tipper bars** -height adjustable
- Seating**
  - tension adjustable straps
  - foam cushion included as standard
  - choice of other cushions available as optional from section 9
- Backrest**
  - tension adjustable straps
  - foam padded upholstery cover included

Optional

TSA117 - FB  
Elevating leg rest

Optional

TSA117 - FG  
Stump support

	overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	backrest height	footrest length	folded width	rear wheel	front castor	unit weight	weight capacity
51cm	min100cm max107cm	106cm	72cm	min40cm max50cm	51cm	51cm	min21cm max28cm	min48cm max53cm	min30cm max39cm	33cm	60cm	18cm	17.4kg	140kg
46cm	min100cm max107cm	106cm	67cm	min40cm max50cm	46cm	51cm	min21cm max28cm	min48cm max53cm	min30cm max39cm	33cm	60cm	18cm	17.4kg	140kg
40cm	min100cm max107cm	106cm	62cm	min40cm max50cm	40cm	51cm	min21cm max28cm	min48cm max53cm	min30cm max39cm	33cm	60cm	18cm	17.4kg	140kg

# PREMIUM PRESCRIPTION WHEELCHAIR - Aluminium



A. **SEAT:** Available in three seat widths. Please select the seat width of your choice.

	Seat width 51cm <input type="checkbox"/>	Seat width 46cm <input type="checkbox"/>	Seat width 40cm <input type="checkbox"/>	\$
---	--	--	--	----

B. **SEAT:** Available in four seat depths. Please select the seat depth of your choice.

	Seat Depth	40cm <input type="checkbox"/>	42.5cm <input type="checkbox"/>	45cm <input type="checkbox"/>	47.5cm <input type="checkbox"/>	50cm <input type="checkbox"/>	\$
---	------------	-------------------------------	---------------------------------	-------------------------------	---------------------------------	-------------------------------	----

C. **CUSHION:** Supplied with a 3cm thick medium density foam cushion as standard. However, there is a variety of other wheelchair cushions available in the Freedom Healthcare range that may better suit your specific needs. Please choose the wheelchair cushion you prefer.

	Seat Cushion, Medium density foam 3cm thick come as standard	\$
	Preferred Seat Cushion <input type="checkbox"/> Product Code.....	\$

D. **SEAT HEIGHT:** Refer to the Seating Diagram Overleaf for available seat heights for this model.

	Floor to Seat Height at Front	48cm <input type="checkbox"/>	49.5cm <input type="checkbox"/>	51cm <input type="checkbox"/>			
	Floor to Seat Height at Rear	42cm <input type="checkbox"/>	44.5cm <input type="checkbox"/>	47cm <input type="checkbox"/>	49.5cm <input type="checkbox"/>	52cm <input type="checkbox"/>	\$

E. **BACKREST:** Height adjustable Backrest. Please fill in your preferred backrest height.

	Backrest Height	48cm <input type="checkbox"/>	50.5cm <input type="checkbox"/>	53cm <input type="checkbox"/>	\$
---	-----------------	-------------------------------	---------------------------------	-------------------------------	----

F. **BACKREST:** Available with an angle adjustable backrest. Please select your preferred backrest angle.

	Backrest Angle	90deg <input type="checkbox"/>	95.6deg <input type="checkbox"/>	101.2deg <input type="checkbox"/>	106.8deg <input type="checkbox"/>	112.4deg <input type="checkbox"/>	\$
---	----------------	--------------------------------	----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	----

G. **ARMREST:** Height adjustable Armrest. Please fill in your preferred Armrest height. Choose any height between the minimum and maximum.

	Armrest Height	20cm <input type="checkbox"/>	21cm <input type="checkbox"/>	22cm <input type="checkbox"/>	23cm <input type="checkbox"/>	24cm <input type="checkbox"/>	25cm <input type="checkbox"/>
		26cm <input type="checkbox"/>	27cm <input type="checkbox"/>				\$

H. **LEGREST:** Comes with swingaway/removable legrests as standard, heel straps included. Also available with elevating legrests or stump support as optional. Please choose your preferred legrest.

	Standard Swingaway/Removable Left Legrest <input type="checkbox"/>	Left Elevating legrest <input type="checkbox"/>	Left Stump Support <input type="checkbox"/>	
	Standard Swingaway/Removable Right Legrest <input type="checkbox"/>	Right Elevating legrest <input type="checkbox"/>	Right Stump Support <input type="checkbox"/>	\$

I. **FOOTPLATES:** Length adjustable footplate. Please fill in your preferred footplate length. Choose any length between the minimum and maximum.

	Footplate Length	Minimum 30cm	Maximum 39cm	.....cm	\$	
	Footplate Angle Position	20° <input type="checkbox"/>	30° <input type="checkbox"/>	40° <input type="checkbox"/>	50° <input type="checkbox"/>	\$

J. **WHEELBASE:** Available with an adjustable wheelbase to aid with personal maneuverability and balance requirements. Please specify your preferred wheelbase (measured from centre of axles)

	Wheelbase	44cm <input type="checkbox"/>	46.5cm <input type="checkbox"/>	49cm <input type="checkbox"/>	51.5cm <input type="checkbox"/>	\$
---	-----------	-------------------------------	---------------------------------	-------------------------------	---------------------------------	----

K. **TYRES:** Available with three tyre options.

	Pneumatic tyres (standard) <input type="checkbox"/>	Airless "Soft Ride" PU tyres <input type="checkbox"/>	High pressure low resistance air tyre <input type="checkbox"/>	\$
---	---	---	--	----

**HANDLE BAR:** Available with height adjustable handle bars. Please fill in your preferred height.

	Handle Bars	Minimum 100cm	Maximum 107cm	.....cm	\$
---	-------------	---------------	---------------	---------	----

N. **LABOUR COSTS:** The above mentioned adjustments may incur a labour cost

	Total labour cost	\$
---	-------------------	----

O. **OPTIONAL EXTRAS:**

	Wheelguards <input type="checkbox"/>	\$
---	--------------------------------------	----

	Lap belt <input type="checkbox"/>	\$
---	-----------------------------------	----

	Wheelchair tray table <input type="checkbox"/>	\$
---	--	----

	Ivy Pole <input type="checkbox"/>	\$
---	-----------------------------------	----

	Oxygen bottle holder <input type="checkbox"/>	\$
---	---	----

	Other <input type="checkbox"/>	\$
---	--------------------------------	----

P. **TOTAL COST:**

	Total cost	\$
---	------------	----

