


Heavy Duty Steel WHEELCHAIR



- Suitable for hospitals, nursing homes and other institutional use.
- Equipped with an array of important features.
- Heavy duty construction.
- Tyres crafted from durable, soft polyurethane.
- No tyre inflation required.
- Designed to fit many optional accessories



Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166
 22 / 170 North Road, Woodridge, QLD 4119 | Ph:1300 460 070
 wecare@mobilityrentals.com.au | www.mobilityrentals.com.au



seat width 51cm

overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	backrest height	footrest length	folded width	rear wheel	front castor	unit weight	weight capacity
min 98cm max 108cm	107cm	69cm	43cm	51cm	51cm	23cm	40cm	min33cm max44cm	32cm	propel 60cm transit 30cm	18cm	19.7kg	130kg

seat width 46cm

overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	backrest height	footrest length	folded width	rear wheel	front castor	unit weight	weight capacity
min 98cm max 108cm	107cm	64cm	43cm	46cm	51cm	23cm	40cm	min33cm max44cm	32cm	propel 60cm transit 30cm	18cm	19.7kg	130kg