


REBOTEC

QUADRO FUN-KIDS

Height adjustable walking frame



- Compression fitting for strong and confident stride without weak weld points
- Wide and stable footprint with 4 interdependent rubber ferrules
- Flexible height with 10 adjustment settings

MPN

Quadro Fun-Kids
156.95

Accessories

Handle Pads (Small)
160.01.01x2



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 297302

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

HEIGHT FROM HANDLE TO CUFF 16 – 21 cm

HEIGHT OF HANDLE 67 – 88 cm

HEIGHT ADJUSTMENTS x10

