

# REBOTEC

## FOX & FIXI

### Anterior walker for children with brakes



Colorful rollator wheeled walker for kids

- Designed as a junior posterior walker specifically for children
- Lightweight & folds easily for storage and travel
- Designs allows full and easy access into the walking area promoting upright posture
- Fitted with a dual brake action, puncture proof tyres and large swivel castors for maximum safety and manoeuvrability

#### MPN

Fox  
186.00.20

Fixi  
186.00.10

#### ACCESSORIES

Back PU Cushion  
280.10.48

Seat, black, wood  
286.65.25

	Fox	Fixi
DEPTH	68cm	68cm
WIDTH	58cm	58cm
WEIGHT	6.3kg	6.2kg
HEIGHT ADJUSTMENT STEPS	5	5
HEIGHT OF HANDLE	70-78cm	61-70cm
MAX BODY WEIGHT	100kg	100kg



Australian Government  
Department of Health  
Therapeutic Goods Administration

ARTG 305299

info@rebotec.com.au  
rebotec.com.au  
218/354 Eastern Valley Way  
Chatswood 2067 NSW

