

# Sturdy Steel WHEELCHAIR



- Suitable for hospitals, nursing homes and other institutional use.
- Equipped with an array of important features.
- Heavy duty construction.
- Tyres crafted from durable, soft polyurethane.
- No tyre inflation required.

**Available with:** - seat width 46cm or 51cm



Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166  
 22 / 170 North Road, Woodridge, QLD 4119 | Ph:1300 460 070  
[wecare@mobilityrentals.com.au](mailto:wecare@mobilityrentals.com.au) | [www.mobilityrentals.com.au](http://www.mobilityrentals.com.au)



**1 YEAR WARRANTY**  
 \* conditions apply

## seat width 51cm

overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	backrest height	footrest length	folded width	rear wheel	front castor	unit weight	weight capacity
90cm	107cm	68cm	40cm	51cm	51cm	21cm	40cm	min30cm max39cm	31cm	propel 60cm transit 30cm	20cm	19kg	150kg

## seat width 46cm

overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	backrest height	footrest length	folded width	rear wheel	front castor	unit weight	weight capacity
90cm	107cm	63cm	40cm	46cm	51cm	21cm	40cm	min30cm max39cm	31cm	propel 60cm transit 30cm	20cm	19kg	150kg