

SUPERLITE 46cm Aluminium Wheelchair




- Quality double braced aluminium frame.
- Easily transportable, super light 12kg Net weight.
- Compact size when folded.
- Tyres crafted from durable, soft polyurethane.
- No tyre inflation required.

Seat width 46cm



Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166
 22 / 170 North Road, Woodridge, QLD 4119 | Ph:1300 460 070
 wecare@mobilityrentals.com.au | www.mobilityrentals.com.au



overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	backrest height	footrest length	folded width	rear wheel	front castor	unit weight	weight capacity
95cm	107cm	64cm	43cm	46cm	50cm	22cm	45cm	min30cm max40cm	29cm	propel 60cm transit 30cm	15cm	12.7kg	110kg