

SUPERLITE 46cm Aluminium Wheelchair






- Quality double braced aluminium frame.
- Easily transportable, super light 12kg Net weight.
- Compact size when folded.
- Tyres crafted from durable, soft polyurethane.
- No tyre inflation required.

Seat width 46cm



1 YEAR WARRANTY
* conditions apply



Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166
22 / 170 North Road, Woodridge, QLD 4119 | Ph:1300 460 070
wecare@mobilityrentals.com.au | www.mobilityrentals.com.au

overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	backrest height	footrest length	folded width	rear wheel	front castor	unit weight	weight capacity
95cm	92cm	64cm	43cm	46cm	50cm	21cm	45cm	min30cm max40cm	29cm	30cm	15cm	11.8kg	110kg