

# REBOTEC

## WALK-ON KIDS WITH ROLLERS

Height adjustable walking frame



Colorful walking frame for kids with 2 front roller wheels.

- Smooth and stylish powder coated aluminium finish
- Non-slip hand-grips for safety and confidence
- Height adjustable walking frame
- Foldable

### MPN

Small  
175.20

Extra Small  
175.10

	Extra Small	Small
DEPTH	40 cm	40 cm
WIDTH	50 cm	50 cm
HEIGHT OF HANDLE	62-71 cm	70-81 cm
WEIGHT	3.3 kg	3.4 kg
HEIGHT ADJUSTMENTS	x5	x5



Australian Government  
Department of Health  
Therapeutic Goods Administration

ARTG 305299

info@rebotec.com.au  
rebotec.com.au  
218/354 Eastern Valley Way  
Chatswood 2067 NSW

