

AIR COMFORT CHAIR / TILT BED

AIR PRESSURE MANAGEMENT



- The seat and backrest are fitted with an air pressure relief system.
- Adjustable seat (multi-adjust).
- Adjustable backrest (multi-adjust).
- Adjustable legrest (multi-adjust).
- Detachable/Adjustable headrest.
- Locking castors.
- Full metal frame.
- Proven pressure relief for long-term seating.
- Varied seating and sleeping positions.
- Replacable cover.
- Easy to clean.
- Easy to repair and maintain.


Pressure relief seating
 This chair offers the comfort of air cushion seating. Air bladders are fitted in strategic places to minimise the development of bed sores (bed ulcers).

- Fabric**
- surecheck 2-way stretch fabric
 - ultra soft and pliable
 - water proof
 - anti-bacterial
 - fire retardant
 - odour repelling



• Patient reclined position and tilt in space function

• Fitted with four pressure relieving air cushions that distribute the patient's weight evenly and assist in preventing bed sores (bed ulcers)



Air pressure management system ensures even distribution of the patient's weight

Incontinence cover available as an option



overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
135cm	102cm	77cm	50cm	47cm	70cm	78cm	45kg	180kg