

# Heavy Duty Walking Tutor

- New aluminium oval tube frame construction.
- Safe working load of 180kg, Australian standards certified.
- Fitted with new “light touch” hand brakes with a lock brake mechanism.
- The frame is widened between the rear wheels for maximum walking space.
- The Gutter arm walker is foldable for transportation.
- The oval tube provides a directional height adjustment function (Press button adjustment).
- The handles bars are multi adjustable for angle and fore/aft function.



overall height	overall length	overall width	armrest length	folded width	castor size	unit weight	weight capacity
min103cm max127cm	88cm	65cm	min30cm max43cm	28cm	20cm	12.2kg	180kg

**W**  
1 YEAR WARRANTY  
\* conditions apply

Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166  
22 / 170 North Road, Woodridge, QLD 4119 | **Ph:1300 460 070**  
wecare@mobilityrentals.com.au | www.mobilityrentals.com.au