

REBOTEC

PLUTO

Reverse posterior child walker, gait trainer



- Rear braking system prevents unintended roll back
- Height and with adjustable handlebars
- Steerable front wheel with solid puncture proof tires
- Support from behind to promote natural upright posture when standing and walking
- Supports therapeutic exercise to learn natural movements and increase natural balance

MPN

Pluto frame, small
187.00.10

Pluto seat, small
287.20.10

Pluto frame, large
187.00.20

Pluto seat, large
287.20.20

Pluto front-wheel steerlock
287.10.42

Harness/sling, small
287.50.01

Harness/sling, large
287.50.02

	Small	Large
LENGTH	79cm	85cm
WIDTH	64cm	67cm
HEIGHT OF HANDLE	45-81cm	55-91cm
WEIGHT	8.0kg	8.5kg
MAX BODY WEIGHT	60kg	80kg



Australian Government
Department of Health
Therapeutic Goods Administration

ARTG 305299

info@rebotec.com.au
rebotec.com.au
218/354 Eastern Valley Way
Chatswood 2067 NSW

