

# Standard Rollator



- Robust construction - with a safe working limit of 130kg, Australian standards certified.
- New "light touch" lock brake system.
- New 20cm flat ridged wheels to minimise wheel wobble and increase grip.
- Easy fold rounded backrest for comfort.
- For indoor and outdoor use.
- Can be folded without removing the storage bag.
- New nylon bag prevents visibility of valuable items.  
e.g. Wallet, Handbag, Keys



Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166  
 22 / 170 North Road, Woodridge, QLD 4119 | Ph:1300 460 070  
 wecare@mobilityrentals.com.au | www.mobilityrentals.com.au



overall height	overall length	overall width	seat depth	seating width	seat height	folded width	castor size	unit weight	weight capacity
min82cm max97cm	69cm	62cm	33cm	44cm	56cm	23cm	20cm	8.4kg	130kg