

Tray trolley with Handbrakes



- Easy to manoeuvre in confined places.
- Lockable brake system.
- Removable trays.
- Easy to use locking handbrakes



Handle height	overall height	overall length	overall width	unit weight	weight capacity
min10cm max25cm	min77cm max92cm	56cm	55cm	8.2kg	20kg



Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166
 22 / 170 North Road, Woodridge, QLD 4119 | **Ph:1300 460 070**
wecare@mobilityrentals.com.au | www.mobilityrentals.com.au