

COMPACT LIFT CHAIR

Large

Proven pressure relief utilising the unique Air Comfort Seating System.

FEATURES

- Total surface air pressure system for long term pressure relief
- Removable seat and back cushions
- Electric lift and recline mechanism
- Front and rear castors
- Wall Saver mechanism
- Back-up battery capability
- Easy to clean, repair and maintain
- Available in CarrEx & Macrosuede
- Custom colours available
- Available in single or twin motor

FABRIC



Anti fungal




Fire resistant



Anti bacterial



Wear resistant



SPECIFICATIONS

Description	Overall Width	Overall Height	Overall Length	Seat Depth	Seat Width	Seat Height	Back Height	Max User Weight
Compact Lift Chair Large AC59020	820mm	1060mm	890mm	510mm	600mm	520mm	740mm	140kg
Compact Lift Chair Large Twin Motor AC59066	820mm	1060mm	890mm	510mm	600mm	520mm	740mm	140kg
Compact Lift Chair Large Mobile AC59101	820mm	1060mm	890mm	510mm	600mm	520mm	740mm	140kg

Improving your life everyday™

kcare.com.au
1800 783 783

KCare
Healthcare Solutions