

COMPACT LIFT CHAIR

Small

Proven pressure relief utilising the unique Air Comfort Seating System.

FEATURES

- Total surface air pressure system for long term pressure relief
- Removable seat and back cushions
- Electric lift and recline mechanism
- Front and rear castors
- Wall Saver mechanism
- Back-up battery capability
- Easy to clean, repair and maintain
- Available in CarrEx & Macrosuede
- Custom colours available
- Available in single or twin motor

FABRIC



Anti fungal



Fire resistant



Anti bacterial



Wear resistant



SPECIFICATIONS

Description	Overall Width	Overall Height	Overall Length	Seat Depth	Seat Width	Seat Height	Back Height	Max User Weight
Compact Lift Chair Small AC59018	750mm	1060mm	820mm	440mm	500mm	480mm	680mm	140kg
Compact Lift Chair Small Twin Motor AC59064	750mm	1060mm	820mm	440mm	500mm	480mm	680mm	140kg
Compact Lift Chair Small Mobile AC59100	750mm	1060mm	820mm	440mm	500mm	480mm	680mm	140kg

Improving your life everyday™

kcare.com.au
1800 783 783

KCare
Healthcare Solutions