

MAXI DELUXE V2

Unparalleled pressure redistribution and adjustable positioning

FEATURES

- Recommended patient height 1620 – 1880mm
- Patented Airgonomic Seating System uses four air cells to immerse and cushion the users
- Electric adjustable back recline, seat title and raising leg rest
- Drop away / adjustable arm rests and back wings
- Carrflex extreme multi-stretch covering with sanitised treatment and enhanced durability
- Large castors make the chair easy to manoeuvre
- One directional local castor to assist with straight line tracking
- Access vehicle approved to Australian Standards
- Custom colours available



SPECIFICATIONS

Description	Overall Length	Overall Height	Overall Width	Seat Width	Floor to Seat Height	Product Weight	Max User Weight
Maxi Deluxe V2 AC59115	1360mm	1385mm	950mm	720mm	640mm	100kg	300kg

Improving your
life everyday™

kcare.com.au
1800 783 783

KCare
Healthcare Solutions