

Aluminium Over Toilet Aid (OTA)

- Constructed from quality aluminium tubing components.
- Painted up to 110 microns. That is 3x thicker than normal coating to prevent scratches.
- New “S” bend design provides maximum space for the patient to fit between the arm rests.
- Seat and the lid are made from a high density plastic to minimise flex.
- Offers a higher safe working load.
- Ultra light weight aluminium frame.
- Fitted with rubber suction tips.
- Low maintenance

SEAT WITH LID



- larger seat and lid
- non flexible
- removable for easy cleaning



- Frame •
- design for maximum seat space and comfort
 - unique “A” frame shape

- Thicker, durable powder coated finish to help protect against scratches and scuffs

- Splash guard included

Ultra light aluminium frame

- Extension legs •
- height adjustable

- Suction rubber tips

overall height	overall length	overall width	seat depth	seating width	seat height	armrest height	unit weight	weight capacity
min65cm max80cm	61cm	58cm	43cm	50cm	min44cm max59cm	21cm	3.65kg	140kg

W
1 YEAR WARRANTY
* conditions apply