


Light Weight Aluminium Shower Chair



- Constructed from aluminium tubing.
- Painted up to 110 microns. That is 3x thicker than normal coating to prevent scratches.
- New “S” bend design provides maximum space for the patient to fit between the arm rests.
- The seat is made from a high density plastic to minimise flex.
- This product offers a higher safe working load.
- Super light aluminium frame.
- Fitted with rubber suction tips.
- Low maintenance



Seat ●
- solid plastic
- perforated seat for drainage

Frame ●
- design for maximum seat space and comfort
- unique “A” frame shape

Extension legs ●
- height adjustable

Ultra light aluminium frame

● Thicker, durable powder coated finish to help protect against scratches and scuffs

● Suction rubber tips
- slip resistant

overall height	overall length	overall width	seat depth	seating width	seat height	armrest height	unit weight	weight capacity
min79cm max94cm	64cm	60cm	41cm	49cm	min44cm max59cm	21cm	3.3kg	140kg

