


Multi-Adjust Rollator



- New height adjustable seating design.
- Designed with an easy height adjustment system.
- New “light touch” lock brake system.
- Primarily suitable for indoor use.
- Can be folded without removing the storage bag.
- New nylon bag prevents visibility of valuable items.
- e.g. Wallet, Hand Bag, Keys

Seat - Height adjustment system

the seat height can be set to one of the 3 settings, floor to seat height: 45cm, 53cm, 57cm



overall height	overall length	overall width	seat depth	seating width	seat height	folded width	castor size	unit weight	weight capacity
min76cm max98cm	min62cm max76cm	57cm	33cm	45cm	min45cm mid53cm max57cm	21cm	15cm	6.7kg	130kg

