

# Ramps & Transfer Boards






## 25mm or 50mm Threshold Ramp




- reduces tripping hazards
- slip resistant recycled rubber
- easy and quick installation
- lower cost than cement or metal plate
- meets access and mobility gradient 1:8



### 25mm Threshold Ramp

overall height	overall length	overall width
		
25mm	1010mm	200mm

### 50mm Threshold Ramp

overall height	overall length	overall width
		
50mm	1040mm	400mm




## TRANSFER BOARDS

### Transfer Boards - Standard - Wood

The transfer board assists in patient transfer

- made from heavy duty bonded timber
- surface is finished with a smooth reinforced resin






overall length	overall width	weight capacity
		
76cm	28cm	114kg

### TRANSFER BOARD - Polymer

The transfer board assists in patient transfer

- Made from heavy duty polymer material
- smooth surface for easy transfer



overall length	overall width	weight capacity
		
71cm	21cm	114kg