


# SHOWER STOOL - ALUMINIUM



- Constructed from aluminium tubing.
- New design provides maximum space for the patient to fit between the arm rests.
- This product offers a higher safe working load.
- Ultra light aluminium frame.
- Fitted with slip resistant suction tips.
- Painted up to 110 microns. That is 3x thicker than normal coating to prevent scratches.
- Ideal for rental or domestic use.
- Low maintenance



1 YEAR WARRANTY  
\* conditions apply

overall height	overall length	overall width	seat depth	seating width	seat height	armrest height	unit weight	weight capacity
min67cm max82cm	45cm	58cm	34cm	48cm	min49cm max64cm	21cm	2.8kg	140kg

Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166  
22 / 170 North Road, Woodridge, QLD 4119 | **Ph:1300 460 070**  
wecare@mobilityrentals.com.au | www.mobilityrentals.com.au