


Underarm Crutches

Replaceable Soft Rubber Armpads
 Replaceable Soft Rubber Hand Grips
 Simple Light Weight Construction



Armpads ●
 - Soft rubber
 - replaceable

Hand grips ●
 - Soft rubber
 - replaceable

Tall Adult

armpit pad height	armpit to handle	unit weight	weight capacity
min134cm max154cm	min40cm max55cm	0.9kg	110kg

Adult

armpit pad height	armpit to handle	unit weight	weight capacity
min114cm max134cm	min29cm max44cm	0.78kg	110kg

Youth

armpit pad height	armpit to handle	unit weight	weight capacity
min94cm max114cm	min25cm max38cm	0.75kg	110kg

W
 1 YEAR WARRANTY
 * conditions apply

Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166
 22 / 170 North Road, Woodridge, QLD 4119 | **Ph:1300 460 070**
 wecare@mobilityrentals.com.au | www.mobilityrentals.com.au