



# Heavy Duty 56cm Wheelchair

- Equipped with an array of features.
- Heavy duty & robust construction.
- Tyres crafted from durable, soft polyurethane.
- Padded Nylon, spill proof upholstery
- Suitable for in-home use, hospitals, nursing homes and other institutional use.

**WEIGHT CAPACITY  
250KG**



**1 YEAR WARRANTY**  
\* conditions apply

Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166  
22 / 170 North Road, Woodridge, QLD 4119 | **Ph:1300 460 070**  
wecare@mobilityrentals.com.au | www.mobilityrentals.com.au

overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	backrest height	footrest length	folded width	rear wheel	front castor	unit weight	weight capacity
min 97cm max 108cm	106cm	75cm	44cm	56cm	52cm	23cm	40cm	min33cm max44cm	32cm	60cm	20cm	23.3kg	250kg