



# Moderna 3 Fuction Bed

- Quality laminated timber finish
- Contemporary decor and style to fit with a modern look.
- Lift up bed rails/side panels
- Provides a non clinical solution for care needs
- 3 Fuction Electrically adjustable Hospital Bed



Adjustable: head/backrest



Adjustable: bed height



Adjustable: leg/knee bend

overall height	base height (w/o mattress)	overall length	overall width	castor size	unit weight	weight capacity
min 74cm max 126cm	min 20cm max 72cm	212cm	106cm	7.5cm	125kg	180kg