

Single Motor Lift Recline Chair

- Style and design combined with form and functionality.
- This chair offers a high level of comfort for the user.
- Suitable for people requiring a comfortable relaxation chair
- Features include lifting the patient into a standing position and vice versa, the chair also reclines into a semi-bed position for comfortable day naps.

Lift Position



Chair lifts into an upright position, making standing for the patient easier

Recline Position



Reclining function converts the chair into a semi-bed position



- Backrest
 - reclines to a semi-bed position
 - easily removable
- The backrest reclines and the legrest elevates simultaneously
- Armrests
 - wide arms
 - soft padding
- Seat
 - plush padded seating
- Legrest
 - the legrest elevates and the backrest reclines simultaneously

Side pockets •
- elasticised

Dual Battery Backup Transformer



2 x 9 Volt batteries
- not included

Hand Controller



overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
108cm	102cm	92cm	52cm	52cm	50cm	62cm	57kg	180kg



1 YEAR WARRANTY
* conditions apply