

Delux Attendant Propelled Mobile Shower Commode

- Polyurethane (PU) seat, seat-infill, back and armrests are water-impermeable,
- Seat is moulded for comfort and warmth and has a front cut-away seat
- Lightweight, rustproof aluminium frame for ease of manoeuvring
- Simple e-clip height adjustment to suit different toilet heights
- Tool-free height and angle adjustable footplates
- Removable footplates swing in/out to facilitate transfers and are easily
- Armrests flip up for side transfers

Backrest ●
- contoured backrest
- Padded for comfort

Swing away Armrest ●



Seat ●
- soft polyurethane padded
- easy to remove
- easy to clean

Frame ●
- Lightweight, rustproof aluminium frame

Leg rest ●
- swing away and removable



Locking Brakes for 4 wheels



pail with lid included ●

● **Castors**
- all lockable
- low maintenance

Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166
22 / 170 North Road, Woodridge, QLD 4119 | **Ph:1300 460 070**
wecare@mobilityrentals.com.au | www.mobilityrentals.com.au

W
1 YEAR WARRANTY
* conditions apply

| overall height | overall length | overall width | seat depth | seating width | seat height | armrest height | footrest length | castor size | unit weight | weight capacity |
|----------------|----------------|---------------|------------|---------------|-------------|-------------------|----------------------|-------------|-------------|-----------------|
| | | | | | | | | | | |
| 97-102cm | 108cm | 56cm | 43cm | 46cm | 51-56cm | min0cm max31cm | min45cm max52.2cm | 10cm | 17kg | 159kg |