

Quality Homestyle 4 Function Profiling Hospital Bed


Mobility[®]
Rentals & Sales

ndis
REGISTERED
PROVIDER

- Integral side rails included as standard for additional safety
- Dewert® motor, handset and controls ensure highest levels of functionality and reliability
- Easy to use 10 button handset features raised buttons for users with limited dexterity or visual impairment
- Attractive wooden headrest and footboard
- Height adjustable between 40 – 80cm to allow easy care / positioning for users
- Dual motor headrest and knee brake adjustment allows both independent movement as well as synchronised movement of the head and foot rests
- Electric backrest, leg rest, knee brake and lower leg elevation are easy to operate using the handset
- Quick action bolt mechanism allows emergency lowering of the headrest
- Four 10cm (4") braked castors allow secure positioning
- Swivelling castors ensure manoeuvrability over different floor surfaces
- Four section mattress support provides multiple profiling positions
- Easy tool free assembling and dismantling
- Locks into a folded position for easy transportation and storage



4 Function Bed
Incl Trendelenburg



Adjustable: bed height



Adjustable: head/backrest



Adjustable: leg/knee bend



HOME STYLE BED SPECIFICATIONS

WIDTH:	102cm
LENGTH:	213cm
HEIGHT MIN:	40cm
HEIGHT MAX:	80cm
HEAVIEST COMPONENT:	20kg
WEIGHT:	98kg
WEIGHT CAPACITY:	200kg

Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166
22 / 170 North Road, Woodridge, QLD 4119 | **Ph:1300 460 070**
wecare@mobilityrentals.com.au | www.mobilityrentals.com.au

W

2 YEAR WARRANTY
* conditions apply