

QUALITY 4 MOTOR RECLINER LIFT CHAIR - QUAD MOTOR



Style and design combined with quality and functionality. Offers a high level of comfort.

Covered with a choice of premium upholstery, including top grade leather.

Compliant to safe working limit: 150 kg.

High grade durable foam and padding used to provide years of postural support, comfort and safety.

Extra 2 motors provide lumbar and neck support adjustment

Patient Lift Position



Lifts patient into a standing position.

Recline Position



Reclining function converts the chair into a semi-flat bed position.



● Head rest adjusts electrically

- Full Electrically adjustable Lumbar Support
Chair reclines into semi-flat bed position

● Armrests
- soft padding

(49 cm seat width)

● Legrest
- elevating legrest
- controlled independently

- 4 function controller
- easy to use
- controls each chair function separately or in unison
- USB connection on controller

Dual Battery Backup Transformer



Illuminated Hand Controller




FABRIC

LEATHER

Grey



Black



overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
105cm	99cm	81cm	51cm	44cm	51cm	64cm	55kg	150kg