

# Quality Mushroom Fabric Dual Motor Lift Recline Chair

- Soft touch Sofia fabric for added style and comfort
- Memory foam seat upholstery provides added comfort and helps distribute pressure for long periods in chair
- Three tier waterfall back design provides extra comfort and support
- Two motors enable the legrest to be elevated / lowered and backrest to be reclined independently to give more flexible support options
- Backrest can be positioned upright while legrest is elevated, for example to watch TV
- It will rise and tilt users to their feet
- Excellent snooze position for supreme comfort and relaxation
- Quiet motors provide a smooth and discreet operation
- Ergonomically contoured seat and back and padded upholstery for added comfort and support
- Two convenient side pockets for easy storage
- Safety key is built into hand controller to prevent accidental movement and entrapment
- In built fail safe battery back-up system which enables the chair to be raised and lowered in the event of a power failure\*
- Comes with a simple to operate four button handset
- Ideal for individuals who have difficulty standing up and sitting down due to restricted mobility e.g. those with MS, Arthritis, Rheumatism, circulation problems
- Available in 3 colour options – Russet, Mushroom and Graphite
- 1 year warranty\*\*
- Maximum weight capacity 150kg



- Backrest
    - independently reclinable
    - reclines to a semi-bed position
    - easily removable
  - 3 stage controller
    - easy to use
    - controls each chair function separately or in unison
  - Seat
    - plush padded seating
  - Armrests
    - wide arms
    - soft padding
  - Legrest
    - elevating legrest
    - controlled independently
- Chair lifts into an upright position, making standing for the patient easier



Reclining function converts the chair into a semi-bed position

## Hand Controller



overall height	overall length	overall width	seat depth	seat width	seat height	armrest height	unit weight	weight capacity
112cm	96cm	92cm	55cm	51cm	51cm	65cm	60kg	150kg

## Battery Backup Power Pack



6 x 9 Volt batteries not included



## Quality Fabric Upholstery

Mobility Rentals & Sales Southside Pty Ltd ABN:37137136166  
22 / 170 North Road, Woodridge, QLD 4119 | Ph:1300 460 070  
wecare@mobilityrentals.com.au | www.mobilityrentals.com.au