

STANDARD DELUXE V2

Market leader in pressure care seating

FEATURES

- Ease of handling for carers and maximum comfort for patients
- Patented Airogonomic Seating System uses four air cells to immerse and cushion the user
- Adjustable back, recline, seat, tilt and raising leg rest
- Drop away / adjustable arms rest and back wings
- Carrex Extreme multi-stretch materials with sanitised treatment and enhanced durability
- Large castors make the chair easy to manoeuvre
- One directional lock castors to assist with straight line tracking
- Access vehicle approved to Australian standards
- Custom colours available



SPECIFICATIONS

Description	Overall Length	Overall Height	Overall Width	Seat Width	Floor to Seat Height	Product Weight	Max User Weight
Standard Deluxe V2 AC59111	1130mm	1260mm	680mm	470mm	570mm	45kg	180kg

Improving your
life everyday™

kcare.com.au
1800 783 783

KCare
Healthcare Solutions