



# SBI

Dimensions = Millimetres



Shower base height  
Standard 120mm - 140mm  
Use accessory SBR100 up to 140mm - 180mm  
Use accessory SBR150 up to 180mm - 230mm  
Use accessory SBR200 up to 230mm - 280mm

Toilet clearance width between frame is 295mm

405

Can be customised  
see accessories

500 (455 wide)


475

120 - 140

## Product features:

- Weight capacity - 250kg
- Rolls over most standard toilets
- Low maintenance
- Arms - removable / lockable / flip-back
- Height adjustable
- Safety mechanisms
- No rust aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Rotates 360 degrees inside shower (subject to space)
- Footrest - fold-up
- Safety belts - lap & chest belts
- Caster Wheels - 125mm x 4 locking
- Commode system
- Tool-free assembly

Bird's eye view



A = 175

B = 105

C = 345

## Accessory items

- Refer to the accessory pages or on our website