

# SB1 SHOWERBUDDY

## For transfer to a shower with tray step

The ShowerBuddy SB1 shower transfer system will turn your existing bathroom into an accessible bathroom allowing you to enjoy showering again without the excessive cost, delays and stress of a full bathroom remodel.

If your shower has a step or upstand preventing roll-in access, our ShowerBuddy SB1 will easily slide you, in a seated position, across the step and into the shower. Once transferred into the shower the bridge and wheeled base can be disconnected and rolled away, allowing the shower door to close for a more enjoyable shower experience. This chair can be rotated when inside the shower allowing you to face the water outlet or for an assistant to help with cleaning.

### Product Features:

- Weight capacity – 250kg
- Rolls over most standard toilets
- Low maintenance
- Arms - removable / lockable / flip-back
- Height adjustable
- Safety mechanisms
- No rust aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Rotates 360 degrees inside shower (subject to space)
- Footrest - fold-up
- Safety belts - lap & chest belts
- Caster Wheels – 125mm x 4 locking
- Commode system
- Tool-free assembly

### Product Dimensions:

- Chair Width - over caster wheels 520mm
- Chair Width – over arms 570mm
- Chair Width – inside arms 475mm - smaller users see accessory pages
- Chair Length - including footrest 800mm (640mm excluding) Backrest Width 460mm
- Backrest Height from Cushion to top 440mm
- Toilet Height Clearance (maximum) 580mm
- Toilet width clearance at narrowest 295mm (centre of track, approx. 200mm above floor)
- Toilet clearance can depend upon the shape of the toilet pan
- Bridge Length (Standard) 405mm or Accessory 200mm or 300mm or customized, see Accessories
- Seat Cushion Dimensions 460mm W x 430mm D
- Seat Height from Floor 550mm to 675mm (Adjustable)
- Seat Cushion Opening 175mm W x 335mm L
- Shower Base unit dimensions 500mm x 455mm
- Box dimensions L 760mm x W 570mm x H 630mm
- Box Gross Weight 29.7kg

