



Etac Molift Mover 180

A mobile lift for all situations

molift[®]
by Etac

A mobile lift for all situations

The Molift Mover 180 is a compact and light lifter, suitable both in the home care and institutional environment.



An everyday favorite

The Molift Mover 180 has a lightweight alloy construction that makes the lifter very easy to maneuver. The large lifting interval (270 - 1630 mm / 10.6" - 64.2") provides a high lift while also lifting from the floor. Molift Mover 180 has also the option of 2- or 4-point suspension.

Ergonomical handling

For improved caregiver ergonomics, the Molift Mover 180 is equipped with Molift's unique handle that provides multiple and flexible grip possibilities. It is also equipped with Molift's environmentally friendly NiMH-battery.

Molift Mover 180, like all Molift Lifter, is also equipped with the unique Molift Automatic Service Notification system that calculates when service is required and gives a warning signal. The integrated software can also be used to retrieve data about lifter usage and lifting patterns using the Molift Service Tool.

Description

Molift Mover 180

Safe Working Load

180 kg / 396 lbs

Material

Aluminium / Steel

Service Software

Included

Total weight

32 kg / 70.5 lbs

Height

1270 mm / 50"

Length

1260 mm / 49.6"

Width chassis/legs outer

680 mm / 27.2"

Width open legs inner

1010 mm / 39.8"

Lifting height

270 - 1630 mm
10.6" - 64.2"

Lifting interval

1360 mm / 53.6"

Leg height

110 mm / 4.3"

Battery

NiMH 14.4 V - 2.2 Ah

Battery charge time

Ca 3 hours

Lifting speed

40 mm / sec
1.6" / sec

Rated performance

50 lifts at: 75 kg and 500 mm /
165 lbs and 19.7"

Wheel size

75 mm / 3"



Features



Unique handle with many grip possibilities

The unique handle frame provides the best natural grip, creating good ergonomics.



Easy to maneuver

The Molift Mover 180 has a total weight of only 32 kg / 70.5 lbs which, together with the lift's design, makes it easy to maneuver.



Integrated service software

Automatic service notification system determines when the lifter requires service and maintenance.



Many options

The Mover 180 can be used with either 2-point or 4-point slingbar. The slingbars are available in five sizes and can be replaced easily without the use of tools.



Easy lifting from chair

Molift Mover 180 lifts easily from chair, wheelchair, toilet and toilet chair.



Easy lifting from bed

Molift Mover 180 lifts easily from bed with Molift Sling solution or the Molift Stretcher.



Molift Mover 180 has a large lifting range (270 - 1630 mm / 10.6" - 64.2") that allows lifting from the floor as well as high lifts.



Molift Mover 180 in combination with Molift Slings provides a comfortable and more spacious and correct position during the entire transfer. For more information, please visit www.molift.com



The Molift Stretcher is especially designed for lateral transfer of clients with suspected spinal injuries and fractures. It is approved and compatible with X-Ray, CT-scan and MRI scanner.

Accessories



Weight Scale for Molift Mover 180

This scale is approved Class III and complies with all technical standards. The scale can also calculate BMI (Body Mass Index) with a simple touch. With the quick-release coupling, the weight scale can be retro-fitted.

Item no: 1840000



Molift Suspension options

Molift Mover 180 can be used with a 2- or 4-point suspension, and there are 5 variants available.

Standard 4-point Medium	440 x 280 x 235 mm (17.3" x 11" x 9.3")	1830011
4-point Large	550 x 410 x 260 mm (21.7" x 16.2" x 10.3")	1830010
2-point Large	590 mm (23.2")	1830001
2-point Medium	490 mm (19.3")	1830002
2-point Small	390 mm (15.4")	1830003



Molift Stretcher

The Molift Stretcher is especially designed for lateral transfer of clients with suspected spinal injuries and fractures.

Molift Stretcher, standard		
excl lifting loops and safety belts		2150103
Molift Stretcher, with MRI-compatibility		
excl lifting loops and safety belts		2150108
Set of lifting loops (4 pcs)		2150357
Set of straps with auto strap collectors (4 pcs)		2150356
Set of safety belts (3 pcs)		2150355