

Action 4NG HD

STABILITY MEETS STYLE



The Action 4NG HD has outstanding functionality with a user weight of 160kg, offering flexibility in a chair to suit almost everyone.



Making Life's
Experiences Possible

Yes, you can.®

Action 4NG HD

A comfortable, foldable, lightweight wheelchair, the Action 4NG is ideally suited for everyday use.

A key feature of this next-generation model is the innovative Invacare hallmark of cross-compatibility; a shared platform with other models in the Action NG family means that users and advisors can easily update or refurbish their wheelchairs and enjoy maximum opportunities for personalised adjustment.



Features



Tech Specs



Seat height
51cm



Max. safe slope
9°



Total weight unloaded
14kg



Seat width
50.5cm, 55.5cm,
60.5cm




Turning radius
85cm



Transport weight
8kg



Seat depth
45-50cm



Total width
seat width + 20cm



Max. user weight
160kg



Backrest height
43cm



Total length
117cm



Crash Test Approved
Fulfills the requirements of ISO 7176-19

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

All rights reserved. All information quoted is believed to be correct at time of publication. Invacare reserves the right to alter product specifications without prior consent. 05/2024



Yes, you can.®