

Safety is a key feature of the Leo, but this does not detract from its stylish and sporty looks. Leo offers users the freedom and confidence to enjoy their essential daily outings and leisure excursions. With high standards of design and safety engineering, the Leo undoubtedly sets a new market standard.



Making Life's Experiences Possible



Leo

The Invacare Leo offers a smooth, safe drive in a scooter that is easy to control and manoeuvre. The Leo intelligent lighting package includes indicators that automatically switch off and a rear brake light for safe driving. In addition to a visual warning for the battery level display, users are also alerted with an audible warning (beep) when the battery needs recharging.

The Leo is a compact mini-scooter that can be easily dismantled for transportation during longer trips or weekend breaks. With a smooth drive and easy manoeuvrability, getting round tight corners in the local shopping centre, museum or leisure complex becomes refreshingly easy. Access to the battery and electronics for maintenance or servicing is much easier with the Leo as only one part needs to be removed.

Features



Users can easily adjust the tiller to suit their needs.



Rear brake lighting warns other people that the scooter is slowing down, thus preventing collisions.



Swivel functionality and sliding rails ensure optimal seating position.



This unique two-step disengaging lever prevents the scooter from accidentally freewheeling.

Tech Specs



Backrest width 470 mm



Seat depth 410 mm



Floor to seat height 460 - 560 mm



Backrest height 475 mm



Driving unit width 590 mm



Total length



Empty weight 59 kg (exc batteries) 83 kg (inc batteries)



Max user weight



Battery capacity



Turning diameter



Max climbable obstacle height 60 mm



Ground clearance 70 mm



Max safe slope



Max driving distance 38 km



Max driving speed 8 km/h



Motor capac 240 W

Colour Options







* Warranty Terms and Conditions apply. Contact Invacare for full details.

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit **www.invacare.co.nz Australia:** Freephone 1800 460 460, email orders@invacare.com.au or visit **www.invacare.com.au**

