


Küschall® Champion

Price List and Order Form



TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____



Address _____ City _____ State _____ Post Code _____ Country _____

Contact Name _____ Comments _____

	Code	Description	Price
<input type="radio"/>	DHA 0031	Küschall Champion	
<input type="radio"/>	DHA 0000	MyKüschall special construction	<i>Complete all information at end of script</i>
SEATING CHOICE			
FRAME			
<input type="radio"/>	DHA 1150	Frame angle 75°, Aluminium	<i>Not with SW 48cm</i>
<input type="radio"/>	DHA 1160	Frame angle 90°, Aluminium	<i>Not with SW 48cm</i>
<input type="radio"/>	DHA 1170	Frame angle 75°, Aluminium Tapered	<i>Not with SW 36-38cm or SK Option</i>
<input type="radio"/>	DHA 1180	Frame angle 90°, Aluminium Tapered	<i>Not with SW 36-38cm or SK Option</i>
FRAME OPTION			
<input type="radio"/>	DHA 1250	SK Option: Folding Front Frame	<i>Not with tapered frames, SW 48cm, LLL < 38cm</i>
SEAT WIDTH (SW)			
<input type="radio"/>	DHA 1320	36cm	<i>Not with tapered frames or fixation strap</i>
<input type="radio"/>	DHA 1325	38cm	<i>Not with tapered frames or fixation strap</i>
<input type="radio"/>	DHA 1330	40cm	
<input type="radio"/>	DHA 1335	42cm	
<input type="radio"/>	DHA 1340	44cm	<i>Not with light seat upholstery</i>
<input type="radio"/>	DHA 1345	46cm	<i>Not with light seat upholstery</i>
<input type="radio"/>	DHA 1350	48cm	<i>Only with tapered frame; Not with light seat upholstery or SK Option</i>
SEAT DEPTH (SD)			
<input type="radio"/>	DHA 1405	34cm	<i>Not with light seat upholstery</i>
<input type="radio"/>	DHA 1415	36cm	<i>Not with light seat upholstery</i>
<input type="radio"/>	DHA 1425	38cm	
<input type="radio"/>	DHA 1435	40cm	
<input type="radio"/>	DHA 1445	42cm	
<input type="radio"/>	DHA 1460	44cm	<i>Not with light seat upholstery</i>
<input type="radio"/>	DHA 1470	46cm	<i>Not with light seat upholstery</i>
<input type="radio"/>	DHA 1480	48cm	<i>Not with light seat upholstery</i>
BACK HEIGHT (BH)			
<input type="radio"/>	DHA 1615	30cm	<i>Not with standard back upholstery</i>
<input type="radio"/>	DHA 1620	31.5cm	<i>Not with standard back upholstery</i>
<input type="radio"/>	DHA 1625	33cm	<i>Not with standard back upholstery</i>
<input type="radio"/>	DHA 1630	34.5cm	
<input type="radio"/>	DHA 1635	36cm	
<input type="radio"/>	DHA 1640	37.5cm	
<input type="radio"/>	DHA 1645	39cm	
<input type="radio"/>	DHA 1650	40.5cm	
<input type="radio"/>	DHA 1655	42cm	
<input type="radio"/>	DHA 1660	43.5cm	
<input type="radio"/>	DHA 1665	45cm	
<input type="radio"/>	DHA 1670	46.5cm	
FRONT SEAT-TO-FLOOR-HEIGHT (FSTF)			
<input type="radio"/>	DHA 1740	45cm	<i>Not with LLL 43-50cm</i>
<input type="radio"/>	DHA 1745	46cm	<i>Not with LLL 44-50cm</i>
<input type="radio"/>	DHA 1750	47cm	<i>Not with LLL 45-50cm</i>
<input type="radio"/>	DHA 1755	48cm	<i>Not with LLL 46-50cm</i>
<input type="radio"/>	DHA 1760	49cm	<i>Not with LLL 47-50cm</i>
<input type="radio"/>	DHA 1765	50cm	<i>Not with LLL 48-50cm</i>
<input type="radio"/>	DHA 1770	51cm	<i>Not with LLL 49-50cm</i>
<input type="radio"/>	DHA 1775	52cm	<i>Not with LLL 50cm</i>
<input type="radio"/>	DHA 1780	53cm	<i>Not with 3" front wheel</i>
<input type="radio"/>	DHA 1785	54cm	<i>Only with 5" front wheels</i>

NOTE: Measurements of the Küschall Champion may not be similar to previous generations - make sure to always take measurements

Küschall® Champion


REAR SEAT-TO-FLOOR-HEIGHT (RSTF)		75°: RSTF=[FSTF-7cm] to [FSTF-0cm]; 90°: RSTF=[FSTF-10cm] to [FSTF-3cm]		
<input type="radio"/>	DHA 1815	39cm	Not with 25" or 26" rear wheels	RSTF / Rear wheel size 22": 39-44 cm 24": 39-46 cm 25": 41-48 cm 26": 43-48 cm
<input type="radio"/>	DHA 1820	40cm	Not with 25" or 26" rear wheels	
<input type="radio"/>	DHA 1825	41cm	Not with 26" rear wheels	
<input type="radio"/>	DHA 1830	42cm	Not with 26" rear wheels	
<input type="radio"/>	DHA 1835	43cm		
<input type="radio"/>	DHA 1840	44cm		
<input type="radio"/>	DHA 1845	45cm	Not with 22" rear wheels	
<input type="radio"/>	DHA 1850	46cm	Not with 22" rear wheels	
<input type="radio"/>	DHA 1855	47cm	Not with 22" or 24" rear wheels	
<input type="radio"/>	DHA 1860	48cm	Not with 22" or 24" rear wheels	
LOWER LEG LENGTH (LLL)		75°: LLL ≤ [FSTF-3cm]; 90°: LLL ≤ [FSTF-4cm]		
<input type="radio"/>	DHA 1924	32cm	Not with SK Option	<p>User's LLL must be measured by taking away the cushion thickness</p>  <p>LLL = LLLpat - SK</p>
<input type="radio"/>	DHA 1926	33cm	Not with SK Option	
<input type="radio"/>	DHA 1928	34cm	Not with SK Option	
<input type="radio"/>	DHA 1930	35cm	Not with SK Option	
<input type="radio"/>	DHA 1932	36cm	Not with SK Option	
<input type="radio"/>	DHA 1934	37cm	Not with SK Option	
<input type="radio"/>	DHA 1936	38cm	SK Option only with seat depth (SD) ≥ 40cm	
<input type="radio"/>	DHA 1938	39cm	SK Option only with seat depth (SD) ≥ 38cm	
<input type="radio"/>	DHA 1940	40cm		
<input type="radio"/>	DHA 1942	41cm		
<input type="radio"/>	DHA 1944	42cm		
<input type="radio"/>	DHA 1946	43cm		
<input type="radio"/>	DHA 1948	44cm		
<input type="radio"/>	DHA 1950	45cm		
<input type="radio"/>	DHA 1952	46cm		
<input type="radio"/>	DHA 1954	47cm		
<input type="radio"/>	DHA 1956	48cm		
<input type="radio"/>	DHA 1958	49cm		
<input type="radio"/>	DHA 1960	50cm		
BACKREST				
<input type="radio"/>	DHA 2140	Angle adjustable, forward folding backrest with soft locking		
BACK UPHOLSTERY				
<input type="radio"/>	DHA 2200	Standard upholstery	Not with BH < 33cm, height-adj rear-set push handles and backrest stabilisation bar	
<input type="radio"/>	DHA 2220	Tension adjustable, black		
<input type="radio"/>	DHA 2230	Tension adjustable, Aero		
PUSH HANDLES				
<input type="radio"/>	DHA 2400	Without push handles	Not with BH < 34.5cm	
<input type="radio"/>	DHA 2410	Push handles fixed, mini		
<input type="radio"/>	DHA 2420	Push handles fixed, long		
<input type="radio"/>	DHA 2430	Height-adjustable rear set	Not with standard upholstery or tubular armrests	
<input type="radio"/>	DHA 2440	Height-adjustable, integrated	Not with backrest stabilisation bar	
<input type="radio"/>	DHA 2450	Push handles, foldable	Not with backrest stabilisation bar	
BACKREST ANGLE (BA)				
<input type="radio"/>	DHA 2510	76°		
<input type="radio"/>	DHA 2530	80.5°		
<input type="radio"/>	DHA 2550	85°		
<input type="radio"/>	DHA 2570	89.5°		
<input type="radio"/>	DHA 2590	94°		
BACKREST STABILISATION BAR				
<input type="radio"/>	DHA 2595	Backrest stabilisation bar, swing away removable	Not with standard upholstery, folding push handles or BH < 40.5cm	
SEATING				
SEAT UPHOLSTERY				
<input type="radio"/>	DHA 2600	Standard seat upholstery		
<input type="radio"/>	DHA 2630	Light seat upholstery	Not with SW and SD > 42cm, or SD 34cm	
SEATING				
SEAT CUSHION				
<input type="radio"/>	DHA 2700	Standard seat cushion, 60mm, double foam with visco top		<input type="radio"/> CLEAR
<input type="radio"/>		Matrix PS Cushion	Not all seat sizes available. Refer to Seating Catalogue	
<input type="radio"/>		Matrix PSP Cushion	Not all seat sizes available. Refer to Seating Catalogue	
<input type="radio"/>		Matrix VI Cushion	Not all seat sizes available. Refer to Seating Catalogue	
<input type="radio"/>		Matrix Libra Cushion	Not all seat sizes available. Refer to Seating Catalogue	
SEAT BAG				
<input type="radio"/>	DHA 2800	Seat bag, click in and removable	Includes fixation clips DHA7010	

Küschall® Champion

FOOTREST			
<input type="radio"/>	DHA 3240	Flip-up footplate, angle adjustable, set to front	
<input type="radio"/>	DHA 3250	Flip-up footplate, angle adjustable, set to back	
CALF STRAP			
<input type="radio"/>	DHA 3320	Mounted on supporter	
<input type="radio"/>	DHA 3330	Mounted on frame	
SIDE GUARDS			
SIDE GUARDS / CLOTHES GUARDS			
<input type="radio"/>	DHA 3420	Carbon clothes guards, flip back	
<input type="radio"/>	DHA 3430	Carbon mud guards, fixed	Not with armrest DHA3490
<input type="radio"/>	DHA 3440	Carbon mud guards, removable	Not with armrest DHA3490
ARMRESTS			
<input type="radio"/>	DHA 3490	Armrest, height adjustable, short armpad (T-arm)	No side guard included, not with mudguards
<input type="radio"/>	DHA 3480	Tubular armrest, height adjustable, removable, comfort padded	
FRAME COLOUR & CUSTOMISATION			
FRAME COLOUR			
<input type="radio"/>	DHA 4000	Space Black	
<input type="radio"/>	DHA 4010	Anthracite Grey	
<input type="radio"/>	DHA 4020	Pure White	Not with Light Grey Bold Styling
<input type="radio"/>	DHA 4030	Carmine Red	
<input type="radio"/>	DHA 4040	Azure Blue	
<input type="radio"/>	DHA 4100	Black Glossy	Extra lead time 10+ days. Tag styling only
<input type="radio"/>	DHA 4110	Ferrite	Extra lead time 10+ days. Tag styling only
<input type="radio"/>	DHA 4120	Slate Metallic	Extra lead time 10+ days. Tag styling only
<input type="radio"/>	DHA 4130	Silver Grey	Extra lead time 10+ days. Tag styling only
<input type="radio"/>	DHA 4140	Royal Blue	Extra lead time 10+ days. Tag styling only
<input type="radio"/>	DHA 4150	Ruby Red	Extra lead time 10+ days. Tag styling only
<input type="radio"/>	DHA 4160	Smoky Grey	Extra lead time 10+ days. Tag styling only
<input type="radio"/>	DHA 4170	Dynamite Orange	Extra lead time 10+ days. Tag styling only
<input type="radio"/>	DHA 4180	Forest Green	Extra lead time 10+ days. Tag styling only
STYLING			
<input type="radio"/>	DHA 4190	Bold	Not with SK Option
<input type="radio"/>	DHA 4195	Tag	
BOLD COLOUR			
<input type="radio"/>	DHA 4200	Black	
<input type="radio"/>	DHA 4210	Dark Grey	
<input type="radio"/>	DHA 4220	Light Grey	
<input type="radio"/>	DHA 4230	Red	
<input type="radio"/>	DHA 4240	Blue	
<input type="radio"/>	DHA 4250	Gold	
TAG COLOUR			
<input type="radio"/>	DHA 4300	Black	
<input type="radio"/>	DHA 4310	Red	
FRONT WHEEL			
FRONT WHEEL			
<i>For configuration please also see tables at end of form</i>			
<input type="radio"/>	DHA 5000	3" Supersport	Not with FSTF 53-54cm
<input type="radio"/>	DHA 5050	4" Everyday wheel, black	Not with FSTF 54cm
<input type="radio"/>	DHA 5100	4" Comfort wheel, black	Not with FSTF 54cm
<input type="radio"/>	DHA 5150	4" Starec, black	Not with FSTF 54cm
<input type="radio"/>	DHA 5200	4" EPIC Softroll, black	Not with FSTF 54cm
<input type="radio"/>	DHA 5250	4" EPIC Softroll, silver	Not with FSTF 54cm
<input type="radio"/>	DHA 5300	5" Everyday wheel, black	
<input type="radio"/>	DHA 5350	5" Comfort wheel, black	
<input type="radio"/>	DHA 5400	5" Starec, black	
<input type="radio"/>	DHA 5450	5" EPIC Softroll, black	
<input type="radio"/>	DHA 5500	5" EPIC Softroll, silver	

Küschall® Champion

REAR WHEELS			
REAR WHEEL POSITION			
<input type="radio"/>	DHA 6100	1st from front	Most active position
<input type="radio"/>	DHA 6110	2nd from front	
<input type="radio"/>	DHA 6120	3rd from front	
<input type="radio"/>	DHA 6130	4th from front	
<input type="radio"/>	DHA 6140	5th from front	
<input type="radio"/>	DHA 6150	6th from front	
<input type="radio"/>	DHA 6160	7th from front	
<input type="radio"/>	DHA 6170	8th from front	
<input type="radio"/>	DHA 6180	9th from front	Least active position
QUAD QUICK RELEASE AXLE			
<input type="radio"/>	DHA 6300	Tetra Clip axle	Not with Spinerger LXX wheels or Alber Add On Kits
REAR WHEEL SIZE			
<input type="radio"/>	DHA 6400	22"	Starec wheels only
<input type="radio"/>	DHA 6410	24"	
<input type="radio"/>	DHA 6420	25"	Not with Starec wheels
<input type="radio"/>	DHA 6430	26"	Starec wheels only
REAR WHEEL TYPE			
<input type="radio"/>	DHA 6510	22", 24", 26" Kuschall Starec, cross spoked	Colourless Starec hub, black rim
<input type="radio"/>	DHA 6520	24", 25" High Performance	Black hub & rim; Not with solid tyres
<input type="radio"/>	DHA 6550	24", 25" Spinerger LXX	Hub, rim & spokes black; Not with solid tyres
REAR WHEEL TYRE			
<input type="radio"/>	DHA 6600	22", 24", 25", 26" Schwalbe Rightrun black/grey	
<input type="radio"/>	DHA 6620	24", 25", 26" Schwalbe Marathon Plus tyre	SmartGuard puncture protected
<input type="radio"/>	DHA 6630	24", 25" Schwalbe ONE	V-Guard puncture resistant
<input type="radio"/>	DHA 6640	22", 24" Profile tyre 7 bar	
<input type="radio"/>	DHA 6650	24" Solid tyre, profile, grey	Not with high performance or Spinerger LXX wheels
HANDRIM			
<input type="radio"/>	DHA 6780	22", 24", 25", 26" Light Alu, grey, dist 3cm	Short tabs with one hole at 3cm
<input type="radio"/>	DHA 6800	22", 24", 25", 26" Light Alu, grey, dist 4cm	Long tabs with two holes at 3cm or 4cm, mounted at 4cm
<input type="radio"/>	DHA 6790	22", 24", 25", 26" Light Alu, black, dist 3cm	Short tabs with one hole at 3cm
<input type="radio"/>	DHA 6700	22", 24", 26" Supergrripp, dist 3cm	Long tabs with two holes at 3cm or 4cm, mounted at 3cm
<input type="radio"/>	DHA 6710	24", 26" Steel polished, dist 3cm	Short tabs with one hole at 3cm
<input type="radio"/>	DHA 6720	24" Titanium, dist 3cm	Short tabs with one hole at 3cm
<input type="radio"/>	DHA 6730	24" Titanium, dist 4cm	Long tabs with two holes at 3cm or 4cm, mounted at 4cm
<input type="radio"/>	DHA 6740	24", 25" Carbolife Gekko, dist 3cm	Long tabs with two holes at 3cm or 4cm, mounted at 3cm
<input type="radio"/>	DHA 6750	24", 25" Carbolife Gekko S, dist 3cm	Long tabs with two holes at 3cm or 4cm, mounted at 3cm
<input type="radio"/>	DHA 6760	24" Carbolife Curve L, dist 3cm	Long tabs with two holes at 3cm or 4cm, mounted at 3cm
<input type="radio"/>	DHA 6770	24" Carbolife Quadro, dist 3cm	Long tabs with two holes at 3cm or 4cm, mounted at 3cm
REAR WHEEL CAMBER			
<input type="radio"/>	DHA 6820	Camber 0°	
<input type="radio"/>	DHA 6830	Camber 1°	Not with Alber Add On Kits
<input type="radio"/>	DHA 6840	Camber 3°	Not with Alber Add On Kits or performance brakes
SPOKE GUARDS			
<input type="radio"/>	DHA 6900	22", 24", 26" Transparent	Not with Spinerger LXX
<input type="radio"/>	DHA 6910	22", 24", 26" Black	Not with Spinerger LXX
ALBER ADD ON KIT			
<input type="radio"/>	DHA 6950	E-Motion & Twion adapter plate	Only with 24" wheels and 0° camber. Includes 2 Alber antitippers
FIXATION CLIPS			
<input type="radio"/>	DHA 7000	Fixation strap	Not with tapered front frame; Not with SW 36-38cm
<input type="radio"/>	DHA 7010	Fixation clip	Not with tapered frame
UNFOLDING HELP			
<input type="radio"/>	DHA 7100	Unfolding help	To open folding mechanism from front
SAFETY			
BRAKES			
<input type="radio"/>	DHA 7200	Pull brake, straight lever	
<input type="radio"/>	DHA 7210	Push brake, bent lever	
<input type="radio"/>	DHA 7220	Push brake, straight lever	
<input type="radio"/>	DHA 7230	Performance brake	Not with 3° camber
<input type="radio"/>	DHA 7240	Active scissor brake	
<input type="radio"/>	DHA 7250	Light Active brake	
BRAKE EXTENSION LEVER			
<input type="radio"/>	DHA 7300	Brake extension lever	Only for straight lever brakes
ANTITIPPER			
<input type="radio"/>	DHA 7500	Antitipper right	
<input type="radio"/>	DHA 7510	Antitipper left	
<input type="radio"/>	DHA 7520	Antitipper both sides	



RSTF / Rear wheel size
 22" : 39-44 cm
 24" : 39-47 cm
 25" : 41-48 cm
 26" : 43-48 cm

CLEAR

CLEAR

Küschall® Champion

SEAT BELT			
<input type="radio"/>	DHA 7800	Seat belt	
PASSIVE LIGHTS (REFLECTORS)			<input type="radio"/> CLEAR
<input type="radio"/>	DHA 7900	Passive lights assembled (mounted on spokes)	<i>Reflectors</i>
<input type="radio"/>	DHA 7910	Passive lights enclosed	<i>Reflectors</i>
ACCESSORIES			
FRAME PROTECTION			
<input type="radio"/>	DHA 8010	Front Frame Protectors, pair	<i>Removable</i>
TIPPER AID <i>Max 2 options per side, eg antitipper, tipper aid, cane holder, transit wheel</i>			
<input type="radio"/>	DHA 8100	Tipper-Aid right	<i>Not on same side as cane holder</i>
CANE HOLDER <i>Max 2 options per side, eg antitipper, tipper aid, cane holder, transit wheel</i>			
<input type="radio"/>	DHA 8200	Cane holder, right	<i>Not on same side as tipper-aid</i>
TRANSIT WHEELS 3"			
<input type="radio"/>	DHA 8300	Transit wheels 3"	<i>Only in combination with anti-tippers both sides</i>
AIR PUMP			
<input type="radio"/>	DHA 8500	Air pump	
TOOL KIT			
<input type="radio"/>	DHA 8550	Tool kit	
WHEEL BAG AND WHEEL KITS			
<input type="radio"/>	DHA 8650	Wheel bag	
<input type="radio"/>	DHA 8680	24" Mountain bike wheels with bag	<i>Mountain bike wheels and tires with alu handrims and quick-release axles</i>
WHEELCHAIR BAGS			
<input type="radio"/>	DHA 8700	Travel bag for Champion	<i>Soft bag, wheels can remain on the wheelchair</i>
<input type="radio"/>	DHA 8750	Travel bag for Champion with SK Option	<i>Rigid bag for folded wheelchair plus luggage</i>

SUPPORT FOR MEASUREMENTS / MYKÜSCHALL REQUESTS

	CUSTOMIZED MEASUREMENTS	
	SW	
	SD	
	BH	
	FSTF	
	RSTF	
	LLL	
	BA	
	CAMBER	
	USER WEIGHT	
USER SIZE		

Küschall® Champion

REAR WHEEL SIZE, BRAKE TYPE, FRAME TYPE, SEAT DEPTH, REAR WHEEL POSITION

	Active Brake 75°/90° frame			Performance brake 75°/90° frame	
	SD 34-40	SD 42	SD 44-48	SD 34-40	SD 42-48
22"	1 to 9	1 to 9	1 to 9	1 to 9	1 to 9
24"	3 to 9	1 to 9	1 to 9	1 to 9	1 to 9
25"	4 to 9	2 to 9	1 to 9	2 to 9	1 to 9
26"	5 to 9	3 to 9	1 to 9	3 to 9	1 to 9

FRAME TYPE, STF DIFFERENCE, FOOTPLATE, SEAT DEPTH, FRONT WHEEL SIZE

FSTF - RSTF	FOOTPLATE TO FRONT						FOOTPLATE TO REAR					
	SD 34-38			SD ≥ 40			SD 34-38			SD ≥ 40		
	3"	4"	5"	3"	4"	5"	3"	4"	5"	3"	4"	5"
75° Frame												
0												
1cm												
2cm												
3cm												
4cm												
5cm												
6cm												
7cm												
90° Frame												
3cm												
4cm												
5cm												
6cm												
7cm												
8cm												
9cm												
10cm												



Local Homologation number, ISO, MDA, TUV logos
 The product includes all standard features
 Maximum user weight: 120kg
 We reserve the right to deliver a + or - tolerance of +/- 10 mm on all measurements
 We reserve the right to make technical changes without notice



NOTE: All specifications and dimensions are approximate.



Invacare Asia-Pacific

Australia:	Freephone 1800 460 460	email orders@invacare.com.au	www.invacare.com.au
New Zealand:	Freephone 0800 468 222	email sales@invacare.co.nz	www.invacare.co.nz
Asia:	Phone +61 2 8839 5330	email exportsales@invacare.com.au	www.invacare.com
Thailand:	Phone +66 0 2108 8555	email csthailandsales@invacare.com	www.invacare.com