

# Introducing the Cobalt 2LiftU Chair Range

**COBALT**  
HEALTH

## 2LIFTU 48cm STANDARD

Code: C48-LIFT



## 2LIFTU 58CM LARGE

Code: C58-LIFT



## 2LIFTU 68CM X-LARGE

Code: C68-LIFT



The 2LiftU is a new type of chair that helps people with mobility challenges get in and out of chairs more easily, while blending in with existing dining settings.

The seat comes in 3 standard widths and is available in a variety of wood finishes and fabric colours to match your décor.

The 2LiftU Chair is made in Europe from high-quality materials designed to last.

An opening at the top of the chair allows for easy handling by users or care givers, and the armrest is designed to provide support when sitting or standing.

- 6 different levels of force
- Reduces the risk of falls
- Improves mobility
- Increases independence
- Wide range of options, including back cushions, anti-decubitus cushions
- 6 standard colours
- 3 wood colours; Beech, Oak, Ash

## CHAIR OPTIONS



2LiftU PLUS lever mechanism allowing raising action to be user controlled available for 48cm wide



3cm high Anti decubitus cushion available in Taupe and Truffle



Height Adjustable backrest cushion available in Taupe, Truffle and Ash

| Chair Width | Code         | Mechanism | SWL   | Seat Height | Seat Depth | Chair Weight | Wood Colours        |
|-------------|--------------|-----------|-------|-------------|------------|--------------|---------------------|
| 48cm        | C48-LIFT-BL  | Standard  | 125kg | 91cm        | 44cm       | 9kg          | Beech<br>Oak<br>Ash |
|             | C48-LIFTP-BL | PLUS      |       |             |            |              |                     |
| 58cm        | C58-LIFT-BL  | Standard  | 250kg |             |            | 10kg         |                     |
| 68cm        | C68-LIFT-BL  | Standard  |       | 11kg        |            |              |                     |

### STANDARD CHAIR FEATURES

- Foam Density: HR30-35kg/m<sup>3</sup>
- 6 different levels of support to help users in and out of the chair, up to 60kg of force
- Nylon, felt or rubber glides under each leg
- Made in the European Union of solid wood
- Fire resistant, anti-bacterial, covers that repel liquids
- Standard chairs with/without arms also available
- Custom colours, fabrics and materials also available in PU, PVC and Fabric Options
- 12 month warranty

### STANDARD FABRIC COLOURS – not all fabrics available with all wood colours



# COBALT

HEALTH

FOR FURTHER INFORMATION OR DEMONSTRATIONS, PLEASE CONTACT US

Victoria

[vic-sales@cobalthhealth.com.au](mailto:vic-sales@cobalthhealth.com.au)

(03) 8521 8444

43 Paramount Blvd, Cranbourne West VIC 3977

New South Wales

[sales@cobalthhealth.com.au](mailto:sales@cobalthhealth.com.au)

(02) 8853 1100

34/9 Salisbury Rd, Castle Hill NSW 2154

PO Box 502 Beecroft NSW 2119

[cobalthhealth.com.au](http://cobalthhealth.com.au)