

Alrick Plexa Bed

Product Code: 7-PLEXA

The Alrick Plexa Bed is a floorline bed designed to provide safe and comfortable sleeping arrangements for users with complex needs. This bed includes full length drop down rails to provide security. It features a locking system, ensuring the rails stay in place while the user is sleeping.



Product Features:

- Australian Made and Owned
- Floorline Model
- Vertical Lift
- Electric Auto-Regression Backrest
- Electric Kneebreak
- Manual Vascular Support
- 250kg SWL
- Full Length Drop Down Rails
- Back Lit Handset
- Square Natural Oak End Boards
- Ironstone Powder Coated Frame



Forward & Reverse Trendelenburg



Chair Position Function



75mm Twin Wheel Castors



Bluetooth Connectivity
handset must be ordered separately

Specifications:

| Size | Single | King Single | Double |
|------------------------------|----------------|----------------|----------------|
| External Dimensions | 1040 x 2150 mm | 1210 x 2150 mm | 1510 x 2150 mm |
| Sleeping Surface | 900 x 2000 mm | 1070 x 2000 mm | 1370 x 2000 mm |
| Low Position | 150 mm | 150 mm | 170 mm |
| High Position | 725 mm | 725 mm | 725 mm |
| Backrest / Knee Angle | 62° / 35° | 62° / 35° | 62° / 35° |
| TR / ATR Angle | 15° / 15° | 15° / 15° | 15° / 15° |
| SWL | 250 kg | 250 kg | 250 kg |



Providing a safe & comfortable sleeping arrangement for users with complex needs.

Optional Features – Added at time of manufacture

- Electric Vascular Support
- Dual Locking Castors
- Bluetooth Handset
- Battery Back Up
- Under Bed Light
- USB Charging Ports
- Custom Powder Coat Colours
- Custom End Board Styles and Colours
- Custom Length Sleeping Surface (2m+)
- Mesh-TEX Sides

Stacked Endboard Colours & Styles

Custom Endboard, styles and colours also available.



Laminate
Square Headboard & Footboard



Upholstered
Plain Headboard & Footboard



Natural Oak
Laminate



Ocean
Upholstered



Ash
Upholstered

Accessories – Added after manufacture

- Mattresses
- Self Help Pole
- Emergency Bed Power Unit
- IV Pole
- EN T-Bar Buffer
- Reading Lamp
- Padded Side Protectors



Custom Jersey
Endboards

