



Etac shower stools

Edge, Easy and Smart shower stools and accessories, designed to assist users and carers with every day routines.

 **etac**[®]

Creating Possibilities

Etac Easy shower stools

Etac Easy are comfortable and safe shower stools that suits most users.

Comfortable seating

The stools have soft contoured seats for the best comfort. The surfaces have an anti-slip pattern, which effectively drains off water so that users always sit safely.

Increased mobility

Easy Swivel allows the user to rotate on the shower stool with little effort and with minimal strain on the body. It is convenient to clean the swivel plate since it can easily be disassembled.



Etac Easy Swivel

Etac Easy



Easy to fit

Practical lock function makes it easy to both fit and remove the legs.



Easy and safe

Height adjustment has built-in safety. The push button locks against the hole when weight is put on the shower stool.



Stability

The legs of the stool are angled outwards to increase stability. This means that the higher seat height, the greater the support base. The soft ferrules adjust perfectly to the floor.



Comfortable and safe seats

The concave seats are comfortable to sit on. The flexibility in the leg fixture entails that the stools can handle any weight displacement which can occur when you move your upper body.

Description	Colour	Width	Depth	Seat Height	Weight	Item no.
Easy	Grey	ø 40 cm 15¾"	38cm 15"	49.5 – 53.5 cm 19½ – 21"	2.8 kg 6.2 lbs	81901010
Easy Swivel	Grey	ø 40 cm 15¾"	38cm 15"	49.5 – 53.5 cm 19½ – 21"	2.8 kg 6.2 lbs	81901050

Construction

Seat: Polypropylene.
Legs: Aluminium & polypropylene



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.