


Foam Mattresses


COBALT
HEALTH



Crash Mat



Memory Foam Mattress



Pressure Care Mattress



Our range of mattresses suits all needs, including crash mats to reduce the risk of injury from falls, as well as convoluted foam and memory foam mattresses for healthcare use.

Product Specs

Crash Mat - Economy

Order code:

Economy : M15

Dimensions: 1800 x 600 x 40 mm

Features

- Water resistant cover
- Reflective strip around edge for visibility at night
- High density foam to provide better falls protection
- No-gap mattress surface when unfolded
- Carry handles

Crash Mat - Deluxe

Order code:

Deluxe - Standard: M06

Dimensions: 1800 x 600 x 50 mm

Deluxe - Extra Large: M09

Dimensions: 1800 x 700 x 60 mm

Features

- Heavy duty water resistant cover
- Anti-slip coating on underside of mat
- Reflective strip around edge for visibility at night
- Concealed zipper allows for cover to be removed for cleaning
- High density foam to provide better falls protection
- No-gap mattress surface when unfolded
- Carry handles

Pressure Care Mattress

Order code:

Single: M11

Dimensions (S): 1980 x 880 x 150 mm

Features

- Firm feel - Multi-layer, high density foam construction
- Top layer of high resilient, diamond pattern foam to provide pressure care for long-term care
- Soft memory-foam heel section helps to protect against pressure ulcers through low skin interface pressure
- Firm and high resilient section along sides to help prevent falls and aid transfers
- Removable and replaceable four-way stretch polyurethane cover
- Water resistant cover
- Fire retardant cover and foam

Memory Foam Mattress

Order codes:

Single: M02

King Single: M02-KS

Dimensions (S): 1980 x 880 x 150 mm

Dimensions (KS): 1980 x 1050 x 150 mm

Features

- Medium firmness
- Multi-layer, high density foam construction
- Top layer of soft memory-foam for comfort and passive pressure relief, protecting against pressure ulcers through low skin interface pressure
- Removable and replaceable four-way stretch polyurethane cover
- Concealed waterfall zip on underside of mattress
- Water resistant cover
- Fire retardant cover and foam

Notice: all products, product specifications and data are subject to change without notice to improve reliability, function, design or otherwise



www.cobalthealth.com.au

sales@cobalthealth.com.au

02 8853 1100

Unit 34, 9 Salisbury Rd, Castle Hill
NSW 2154 Australia